



www.mortara.com

REF 9515-151-51-ENG Rev C1

X-Scribe

DATA EXCHANGE INTERFACES

ADMINISTRATOR MANUAL

Manufactured by Mortara Instrument, Inc., Milwaukee, Wisconsin U.S.A.



CAUTION: Federal law restricts this device to sale by or on the order of a physician.



Copyright © 2015
by Mortara Instrument, Inc.
7865 N. 86th Street
Milwaukee, Wisconsin 53224

This document contains confidential information that belongs to Mortara Instrument, Inc. No part of this document may be transmitted, reproduced, used, or disclosed outside of the receiving organization without the express written consent of Mortara Instrument, Inc. Mortara is a registered trademark of Mortara Instrument, Inc. X-Scribe is a trademark of Mortara Instrument, Inc. DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information. Microsoft and Windows are trademarks or registered trademarks of Microsoft Corporation. Applies to X-Scribe V3.30 and HX-Gate V1.10

TECHNICAL SUPPORT AND SERVICE

Headquarters

Mortara Instrument, Inc.

7865 North 86th Street
Milwaukee, WI 53224
U.S.A.

Tel: 414.354.1600
Tel: 800.231.7437
Fax: 414.354.4760
Internet: <http://www.mortara.com>

European Union Representative

Mortara Instrument Europe, s.r.l.

(European Headquarters)
Via Cimarosa 103/105
40033 Casalecchio di Reno (BO)
Italy
Tel: +39.051.298.7811
Fax: +39.051.613.3582

Service/Technical Support Group

Mortara Instrument, Inc.

7865 North 86th Street
Milwaukee, WI 53224
U.S.A.
Tel: 414.354.1600
Service: 888.MORTARA
(888.667.8272)
Fax: 414.354.4760
E-mail: techsupport@mortara.com

Sales Support/ Supplies & Accessories

Mortara Instrument, Inc.

7865 North 86th Street
Milwaukee, WI 53224
U.S.A.

Tel: 414.354.1600
Fax: 414.354.4760

Hospital Customers: orders.us@mortara.com
Physician Practice: orderspc.us@mortara.com
U.S. Distribution: orderspc.us@mortara.com

Mortara Instrument Germany

Bonifaciusring 15
45309 Essen
Germany
Tel: +49.201.18 55 69 70
Fax: +49.201.18 55 69 77

Mortara Instrument Netherlands

Postbus 324
5680 AH Best
Industrieweg 160b
5683 CG Best
Netherlands
Tel: +31.499.377310
Fax: +31.499.377908

Mortara Instrument Australia

PO Box 7568
Baulkham Hills NSW 2153
Unit 28, 9 Hoyle Avenue
Castle Hill NSW 2154
Australia
Tel: +61 2 8070 9303
Fax: +61 2 9899 9478

Mortara Dolby UK Ltd.

Units 11 & 12, Scion House
Stirling University Innovation Park
Stirling FK9 4NF
Scotland
Tel: +44.1786.444980
Fax: +44.1786.446630

NOTICES

Manufacturer's Responsibility

Mortara Instrument, Inc. is responsible for the effects on safety and performance only if:

- Assembly operations, extensions, readjustments, modifications, or repairs are carried out only by persons authorized by Mortara Instrument, Inc.
- The electrical installation of the relevant room complies with the requirements of appropriate regulations, and
- The device is used in accordance with the instructions for use.

Responsibility of the Customer

The user of this device is responsible for ensuring the implementation of a satisfactory maintenance schedule. Failure to do so may cause undue failure and possible health hazards.

Equipment Identification

Mortara Instrument, Inc. equipment is identified by a serial and reference number on the back of the device. Care should be taken so that these numbers are not defaced.

Copyright and Trademark Notices

This document contains information that is protected by copyright. All rights are reserved. No part of this document may be photocopied, reproduced, or translated to another language without prior written consent of Mortara Instrument, Inc.

Other Important Information

The information in this document is subject to change without notice.

Mortara Instrument, Inc. makes no warranty of any kind with regard to this material including, but not limited to, implied warranties of merchantability and fitness for a particular purpose. Mortara Instrument, Inc. assumes no responsibility for any errors or omissions that may appear in this document. Mortara Instrument, Inc. makes no commitment to update or to keep current the information contained in this document.

WARRANTY INFORMATION

Your Mortara Warranty

MORTARA INSTRUMENT, INC. (hereafter referred to as “Mortara”) warrants that components within Mortara products (hereafter referred to as “Product/s”) will be free from defects in workmanship and materials for the number of years specified on documentation accompanying the product, or previously agreed to by the purchaser and Mortara, or if not otherwise noted, for a period of twelve (12) months from the date of shipment.

Consumable, disposable or single use products such as, but not limited to, PAPER or ELECTRODES are warranted to be free from defects in workmanship and materials for a period of 90 days from the date of shipment or the date of first use, whichever is sooner.

Reusable product such as, but not limited to, BATTERIES, BLOOD PRESSURE CUFFS, BLOOD PRESSURE HOSES, TRANSDUCER CABLES, Y-CABLES, PATIENT CABLES, LEAD WIRES, MAGNETIC STORAGE MEDIUMS, CARRY CASES or MOUNTS, are warranted to be free from defects in workmanship and materials for a period of 90 days. This warranty does not apply to damage to the Product/s caused by any or all of the following circumstances or conditions:

- a) Freight damage;
- b) Parts and/or accessories of the Product/s not obtained from or approved by Mortara;
- c) Misapplication, misuse, abuse, and/or failure to follow the Product/s instruction sheets and/or information guides;
- d) Accident; a disaster affecting the Product/s;
- e) Alterations and/or modifications to the Product/s not authorized by Mortara;
- f) Other events outside of Mortara’s reasonable control or not arising under normal operating conditions.

THE REMEDY UNDER THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT WITHOUT CHARGE FOR LABOR OR MATERIALS, OR ANY PRODUCT/S FOUND UPON EXAMINATION BY MORTARA TO HAVE BEEN DEFECTIVE. This remedy shall be conditioned upon receipt of notice by Mortara of any alleged defects promptly after discovery thereof within the warranty period. Mortara’s obligations under the foregoing warranty will further be conditioned upon the assumption by the purchaser of the Product/s (i) of all carrier charges with respect to any Product/s returned to Mortara’s principal place or any other place as specifically designated by Mortara or an authorized distributor or representative of Mortara, and (ii) all risk of loss in transit. It is expressly agreed that the liability of Mortara is limited and that Mortara does not function as an insurer. A purchaser of a Product/s, by its acceptance and purchase thereof, acknowledges and agrees that Mortara is not liable for loss, harm, or damage due directly or indirectly to an occurrence or consequence therefrom relating to the Product/s. If Mortara should be found liable to anyone under any theory (except the expressed warranty set forth herein) for loss, harm, or damage, the liability of Mortara shall be limited to the lesser of the actual loss, harm, or damage, or the original purchase price of the Product/s when sold.

EXCEPT AS SET FORTH HEREIN WITH RESPECT TO REIMBURSEMENT OF LABOR CHARGES, A PURCHASER'S SOLE EXCLUSIVE REMEDY AGAINST MORTARA FOR CLAIMS RELATING TO THE PRODUCT/S FOR ANY AND ALL LOSSES AND DAMAGES RESULTING FROM ANY CAUSE SHALL BE THE REPAIR OR REPLACEMENT OF DEFECTIVE PRODUCT/S TO THE EXTENT THAT THE DEFECT IS NOTICED AND MORTARA IS NOTIFIED WITHIN THE WARRANTY PERIOD. IN NO EVENT, INCLUDING THE CLAIM FOR NEGLIGENCE, SHALL MORTARA BE LIABLE FOR INCIDENTAL, SPECIAL, OR CONSEQUENTIAL DAMAGES, OR FOR ANY OTHER LOSS, DAMAGE, OR EXPENSE OF ANY KIND, INCLUDING LOSS OF PROFITS, WHETHER UNDER TORT, NEGLIGENCE OR STRICT LIABILITY THEORIES OF LAW, OR OTHERWISE. THIS WARRANTY IS EXPRESSLY IN LIEU OF ANY OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND THE WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE.

USER SAFETY INFORMATION



Warning:

Means there is the possibility of personal injury to you or others.



Caution:

Means there is the possibility of damage to the device.

Note:

Provides information to further assist in the use of the device.



Warning(s)

- Reference X-Scribe™ user manual for warnings relating to the stress exercise system.
- Possible malfunction risks could be associated when installing 3rd party software. Mortara Instrument, Inc. cannot verify the compatibility of all possible hardware/software combinations.
- Data exchange interfaces are used to exchange data with 3rd party ECG management systems. It is not possible to assure complete compatibility with all possible 3rd party ECG management systems and configurations. It is recommended to contact the 3rd party vendor to ensure device has been verified as compatible with a particular installation of their system.



Caution(s)

- Reference X-Scribe user manual for cautions relating to the stress exercise system.



Note(s)

- Reference X-Scribe user manual for notes relating to the stress exercise system.

EQUIPMENT SYMBOLS

Symbol Delineation



Attention, consult accompanying documents



Do not dispose as unsorted municipal waste. Per European Union Directive 2002/96, requires separate handling for waste disposal according to national requirements



Indicates compliance to applicable European Union directives

TABLE OF CONTENTS

INTRODUCTION SECTION 1

Purpose	1
Architecture	1
HL7	3
DICOM	3
DICOM Storage Events	3
DICOM MPPS Events	3

CONFIGURATION OF INTERFACES SECTION 2

Filenames	5
XML Orders	6
Configuration of HL7	6
Configuration of DICOM	6

INSTALLATION SECTION 3

Install HXGate	9
Open Ports in Firewall	9

TROUBLESHOOTING SECTION 4

Log Files	11
DICOM Echo	11

MORTARA STRESS REPORT XML APPENDIX A

XML Description	13
-----------------------	----

STRESS REQUEST (ORDER) XML APPENDIX B

Stress Request (Order) XML Description	25
--	----

ALTERNATE 1 XML APPENDIX C

Alternate 1 XML Description	29
-----------------------------------	----

DICOM ERROR MESSAGES APPENDIX D

DICOM Error Codes	33
-------------------------	----

Purpose

The X-Scribe stress system offers several ways to exchange information with 3rd party ECG management systems. This manual covers the technical aspects of installing and configuring the various interface options.

- Stress report as PDF file
- Stress report as TIFF file
- Stress report statistics as Mortara XML file
- Stress report statistics as Alternate XML file
- Stress report statistics sent in HL7 message
- Stress report as PDF encapsulated in HL7 message
- Stress report as DICOM®-encapsulated PDF storage
- Stress orders as DICOM Modality Worklist (MWL) request
- Stress orders received in HL7 messages
- Stress orders receives as XML files
- Stress report DICOM storage commitment
- Stress study status reported as DICOM Modality Performed Procedure Step (MPPS)


***NOTE:** This manual may contain screen shots. Any screen shots are provided for reference only and are not intended to convey actual operating techniques. Consult the actual screen in the host language for specific wording.*

Architecture

Figure 1-1 shows the overall architecture of the X-Scribe data exchange interfaces. X-Scribe can be configured to export PDF, TIFF, and XML files into folders monitored by other systems for import. The exported files can also be used by Mortara's HL7 Gateway for generating HL7 results messages. The Microsoft® Windows™ service "HXGate for X-Scribe" is used to generate DICOM messages; the Windows service "HXGateWaitForEvent" listens for confirmations from DICOM storage commitment requests.

The HL7 Gateway is centralized per enterprise and typically runs on its own server; however, the Windows' HXGate and HXGateWaitForEvent services are installed on each X-Scribe system.

Figure 1-1



HL7

HL7 implementations vary widely between products, versions, and installations. As such, Mortara offers its HL7 Gateway which includes a commercial HL7 messaging broker that is configured for each installation. Mortara representatives work directly with the site's IT representatives to configure and test the interface.

The HL7 Gateway uses the Mortara XML files to generate results messages. These messages can include discrete data found in the XML file as well as the displayable PDF or TIFF report.

DICOM

HXGate is used to communicate with 3rd party systems that use the DICOM protocol. It supports Modality Worklist (MWL) queries, PDF storage, storage commitment, and Modality Performed Procedure Step (MPPS) status messages.

DICOM Storage Events

DICOM-encapsulated PDF messages are sent whenever a report is exported from the archive, or when the report is saved and "Export on Save" option is enabled. If Storage Commitment is enabled, a commitment request will be sent at the same time.

DICOM MPPS Events

Events in X-Scribe cause MPPS messages to be sent if MPPS is enabled in HXGate:

- An "In Progress" message is sent when a procedure is started. X-Scribe considers the procedure started when a new exam is selected and the "OK" button is clicked after entering the patient demographics.
- A "Discontinued" message is sent when a procedure is stopped before entering into the exercise phase. X-Scribe considers the procedure discontinued when the "Abort" button is clicked before "Start Exercise" is selected.
- A "Completed" message is sent when a procedure is completed. X-Scribe considers a procedure completed when the "End" button is clicked in the recovery phase.

CONFIGURATION OF INTERFACES

SECTION 2

X-Scribe can export several kinds of files. The file types and locations are configured in the Miscellaneous tab of the User Profile Settings.

Option	Description
Export Directory	The full path where the files will be exported. The default is C:\XS2\REPORTS.
Request Directory	The full path where order XML files are placed. The default is C:\XS2\Requests.
Sub Directories	Check if exported files should be placed into a separate sub directory per test. The sub directories will have the following format: X<PatientID>[<ExamDate><ExamTime>]<FirstName>_<MI>_<LastName>.
Export On Save	Check if files should be automatically exported when the report is saved.
Export Formats	Description
PDF	Check if a PDF of the report should be exported.
Mortara XML	Check if Mortara XML should be exported.
Alternate XML	Check is Alternate XML should be exported.
TIFF	Check if a TIFF of the report should be exported.
HL7	Check if orders will be coming from an HL7 interface or from another external source of XML-format orders. If unchecked, the clinician can pre-enter demographics and orders using the X-Scribe interface.
DICOM	Check if a DICOM interface is being used.

Filenames

Configure filename formats in the “XS2.ini” file. The file is located in the folder “C:\Mortara Instrument Inc\XScribeII\System”. Add a line specifying “ExportFileNameFormat=” in the [Behavior] section. An appropriate file extension will be added automatically, e.g., .xml or .pdf. **EXAMPLE:**

```
[Behavior]
ExportFileNameFormat=<PatientID>[<TestDateTime>]<PatientFirstName>_<PatientLastName>
```

The following variables can be used to specify the filename format:

Variable	Description
<PatientID>	Patient's ID.
<PatientLastName>	Patient's last name.
<PatientFirstName>	Patient's first name.
<PatientMiddleName>	Patient's middle name.
<PatientSex>	Patient's gender: M = male F = female U = Unknown
<TestDateTime>	Start date/time of the test in HL7 yyyyMMddHHmmss format.
<TestDate>	Start date of the recording in HL7 yyyyMMdd format.
<ReportDateTime>	The current (local) date/time when the report is exported, in HL7 yyyyMMddHHmmss format.
<ReportDate>	The current (local) date when the report is exported, in HL7 yyyyMMdd format.
<PatientDOB>	Patient's date of birth, in HL7 yyyyMMdd format.

X-Scribe will not replace incompatible characters used as part of the filename format. For example, if the Patient ID field is used and it contains a / (slash), the resulting filename will be incompatible with the Windows file system and the file export will fail.

X-Scribe overwrites files with the same name. If exported files will remain in the destination folder, verify the filename formats include variables to ensure each filename is unique. For example, use <TestDateTime> so multiple tests for the same patient are named uniquely. Or, if each revision (export) of the same test will be kept, use <ReportDateTime> so each revision is named uniquely.

XML Orders

When X-Scribe is configured for HL7 mode, orders in XML files can be dropped into its *Request Directory* from the Mortara HL7 Gateway, or from other systems that drop the XML files directly. The XML files must be in the format described in Appendix B. Filenames do not matter; however, filenames starting with the reserved words “Dicom” and “Request” will be ignored. File extensions must be .xml. Orders in the Request Directory can be seen when X-Scribe displays the Worklist.

X-Scribe will delete the corresponding XML file when the exam is completed. If the exam is aborted before entering the exercise phase, the XML file will remain in the directory. If an order is canceled before it is used, the corresponding XML file must be removed by the source system generating the XML files. The source system is also responsible for replacing XML files with updated versions whenever an order is updated (e.g., patient demographics change).

Configuration of HL7

HL7 implementations vary widely between products, versions, and installations. As such, Mortara offers its Mortara HL7 Gateway which includes a commercial HL7 messaging broker that is configured for each installation. Mortara representatives work directly with the site’s IT representatives to configure and test the interface.

If HL7 order messages will be received, check the HL7 option in the Miscellaneous tab of the User Profile Settings.

Configuration of DICOM

HXGate is used to communicate with 3rd party systems that use the DICOM protocol. All DICOM communications settings and Modality Worklist (MWL) filters are configured by running HXConfig.exe. This configuration utility is found in the HXGate folder where the service was installed, typically at:

C:\Mortara Instrument Inc\HXGate\HXConfig.exe.

The left side of the configuration table shows settings for the four DICOM services: C-STORE, C-FIND/MWL, MPPS, and Storage Commitment. It is assumed that C-STORE and Storage Commitment SCPs share the same host and port.

Setting	Description
C-STORE	
Host	The network name or TCP/IP address of the Storage SCP server. This is the same SCP for Storage Commitment requests.
Port	The TCP/IP port number used by the Storage SCP.
AE Title	Application Entity Title of the Storage SCP.
Retry (days)	When the Storage SCP is unavailable, the number of days to continue trying to reach the Storage SCP before giving up.

Setting	Description
New Series Instance UID	Check if each version of a report should be stored in a new DICOM Series. If unchecked, each version is stored in the same Series.
MWL (Modality Worklist)	
Host	The network name or TCP/IP address of the MWL SCP server.
Port	The TCP/IP port number used by the MWL SCP.
AE Title	Application Entity Title of the MWL SCP.
MPPS (Modality Performed Procedure Step)	
Enable MPPS	Check to enable MPPS messages.
Host	The network name or TCP/IP address of the MPPS SCP server.
Port	The TCP/IP port number used by the MPPS SCP.
AE Title	Application Entity Title of the MPPS SCP.
Storage Commitment	
Enable Storage Commitment	Check to enable Storage Commitment messages. Messages will be sent to the C-STORE SCP after report is stored.
Receive Event Port	TCP/IP port where HXGateWaitForEvent Windows service listens for results messages.
Wait For Response (days)	Number of days to wait for a result message before giving up.

Station Settings		Description
SCU AE Title		The SCU AE Title used by HXGate in all DICOM transactions.
Institution Name	(0008,0080)	This is the name of the institution and/or department where the Holter exam was performed. This will be stored in the DICOM object.
Station Name	(0008,1010)	DICOM Station Name assigned to this X-Scribe system. This will be stored in the DICOM object.
Modality	(0008,0060)	Typically "ECG".
Scheduled Station Name	(0040,0010)	
Scheduled Procedure Step Location	(0040,0011)	
Current Patient Location	(0038,0300)	
Requested Procedure Location	(0040,1005)	
Scheduled Procedure Step ID	(0040,0009)	
Requested Procedure ID	(0040,1001)	
Accession Number	(0008,0050)	
Scheduled Station AE Title	(0040,0002)	
User Tag		This can be used to define one additional DICOM tag for filtering the MWL.
User Tag Value		
Requested Procedure Description List	(0032,1060)	List of procedure descriptions accepted. These values are not sent to the MWL SCP for filtering, but are used by HXGate to filter the results from the MWL SCP.
Use Institution	(0008,0080)	When checked, the Institution Name will be used to filter the Modality Worklist query.

Installation of HXGate includes installation of the HXGate and HXGateWaitForEvent Windows services. This is necessary to support DICOM. HXGate must be installed on every X-Scribe stress system and Review Station when DICOM is used. If PDF and XML file export, or HL7 will be used without DICOM, HXGate is not necessary.

Install HXGate

1. Install X-Scribe first. If X-Scribe is not installed, do not install HXGate.
2. Insert the “HXGate for X-Scribe” installation CD.
3. If the installer does not automatically launch, use Windows Explorer to view the files on the CD and launch **HXGate (X-Scribe).msi**
4. Follow the prompts to complete the installation.
5. When possible, keep the default folder: C:\Mortara Instrument Inc\HXGate\
6. Reboot Windows after installation is complete.

In X-Scribe, check the DICOM box in the Miscellaneous tab under User Profiles in Settings.

Open Ports in Firewall

If the system where HXGate was installed has firewall software, it must be configured to allow connections to the computer on the TCP/IP port it listens on. The default port for HXGateWaitForEvent is 104, but it can be changed in the HXGate configuration dialog.

Log Files

HXGate creates a log file for each day of the month. Log files are kept for one month and are overwritten on the same day of the next month. The log files are saved in the “logfiles” folder where HXGate is installed, typically C:\Mortara Instrument Inc\HXGate\LogFiles. Examine the log files when troubleshooting communication problems.

HXGate uses a DICOM library from Merge and has its own log. Modify the file C:\Mortara Instrument Inc\HXGate\bin\Merge.ini to turn logging on. Remove the comment symbol, #, in front of the last lines of the file to turn on the various levels of logging. There are comments in the file to indicate what is logged at each level. Messages will be logged into the file C:\Mortara Instrument Inc\HXGate\bin\merge.log. When changes are made to the merge.ini file, the HXGate service must be restarted for the changes to take effect.

The Merge library uses an environment variable called MERGE_INI. If there is some question where Merge is looking for its ini and DLL files, check this variable. If it is not set, Merge will use a default search path (e.g., current working folder, system32, etc.).

DICOM Echo

A DICOM Echo should be performed whenever you need to test the connection to the DICOM SCPs. Open a cmd window, or use Windows Explorer and go to the HXGate installation folder (default is C:\Mortara Instrument Inc\HXGate). There are 3 BAT files starting with “Echo”, one each for CSTORE, MPPS, and MWL. Execute the appropriate file. A command window will display the result of the DICOM Echo.

MORTARA STRESS REPORT XML

XSD file name: **XScribeOutput3_30.xsd**

File location: on the X-Scribe installation CD, in the **misc** folder.

NOTE: When ST levels are reported in millimeters, the gain is assumed to be 10 mm/mV.

XML Description

XML Tag	Description
/StressTest	
./Manufacturer	Name of manufacturer of the system that produced the report. Always "Mortara Instrument, Inc."
./Version	Name and version of the system that produced the report.
./PDF_Path	Full path and filename of the exported report in PDF format.
./TIFF_Path	Full path and filename of the exported report in TIFF format.
/StressTest/PatientDemographics	
./LastName/Value	Patient's last name.
./LastName/Label	
./FirstName/Value	Patient's first name.
./FirstName/Label	
./MiddleName/Value	Patient's middle name.
./MiddleName/Label	
./ID/Value	Patient's primary medical record number.
./ID/Label	
./SecondaryID/Value	Patient's alternate ID. Site-specific use.
./SecondaryID/Label	
./DOB/Value	Patient's date of birth in format displayed to the user.
./DOB/Label	
./DobEx/Value	Patient's date of birth in XML format, yyyy-MM-dd.
./DobEx/Label	
./Age/Value	Patient's age at time of exam.
./Age/Units	Always years .
./Age/Label	
./TargetHR/Value	Target heart rate to be reached in this exam.
./TargetHR/Units	Always BPM (beats per minute).
./TargetHR/Label	
./Gender/Value	Patient's gender. Values can be: <ul style="list-style-type: none">• Male• Female• Unknown
./Gender/Label	

XML Tag	Description
./Race/Value	Patient's race, free text or selected from the customizable pick list configured in RaceList.txt . Factory default choices in English are: <ul style="list-style-type: none"> • Caucasian • Afro American • Black • Asian • Hispanic • Indian • Other
./Race/Label	
./Height/Value	Patient's height at time of exam.
./Height/Units	<ul style="list-style-type: none"> • in = inches • cm = centimeters
./Height/Label	
./Weight/Value	Patient's weight at time of exam.
./Weight/Units	<ul style="list-style-type: none"> • lbs = pounds • kg = kilograms
./Weight/Label	
./Address/Value	Patient's home address. House number and street.
./Address/Label	
./City/Value	Patient's home city. Free text or selected from the customizable pick list configured in CityList.txt .
./City/Label	
./State/Value	Patient's home state.
./State/Label	
./PostalCode/Value	Patient's home postal (zip) code.
./PostalCode/Label	
./Country/Value	Patient's home country.
./Country/Label	
./HomePhone/Value	Patient's home phone number.
./HomePhone/Label	
./WorkPhone/Value	Patient's work phone number.
./WorkPhone/Label	
./ReferringPhysician/Value	Full name of referring physician. Free text or selected from the customizable pick list configured in ReferringPhysicianList.txt .
./ReferringPhysician/Label	
./AttendingPhysician/Value	Full name of the attending physician. Free text or selected from the customizable pick list configured in AttendingPhyList.txt .
./AttendingPhysicianLabel	
./Smoker/Value	Values can be: <ul style="list-style-type: none"> • Yes • No • Unknown
./Smoker/Label	

XML Tag	Description
./Diabetic/Value	Values can be: <ul style="list-style-type: none"> • Yes • No • Unknown
./Diabetic/Label	
./HistoryOfMI/Value	Values can be: <ul style="list-style-type: none"> • Yes • No • Unknown
./HistoryOfMI/Label	
./FamilyHistory/Value	Values can be: <ul style="list-style-type: none"> • Yes • No • Unknown
./FamilyHistory/Label	
./PriorCABG/Value	Values can be: <ul style="list-style-type: none"> • Yes • No • Unknown
./PriorCABG/Label	
./PriorCath/Value	Values can be: <ul style="list-style-type: none"> • Yes • No • Unknown
./PriorCath/Label	
./Angina/Value	Values can be: <ul style="list-style-type: none"> • Atypical • Typical • None
./Angina/Label	
./Indications/Line	There is one Line per indication. Free text or selected from the customizable pick list configured in IndicationsList.txt .
./Indications/Label	
./Medications/Line	There is one Line per medication. Free text or selected from the customizable pick list configured in MedicationsList.txt .
./Medications/Label	
./Notes/Line	There is one Line per note. Free text or selected from the customizable pick list configured in NotesList.txt .
./Notes/Label	

XML Tag	Description
./MessageID/Value	Values carried over from the XML requests.
./MessageID/Label	
./OrderNumber/Value	
./OrderNumber/Label	
./BillingCode1/Value	
./BillingCode1/Label	
./BillingCode2/Value	
./BillingCode2/Label	
./BillingCode3/Value	
./BillingCode3/Label	
./ExpansionField1/Value	
./ExpansionField1/Label	
./ExpansionField2/Value	
./ExpansionField2/Label	
./ExpansionField3/Value	
./ExpansionField3/Label	
./ExpansionField4/Value	
./ExpansionField4/Label	
/StressTest/TestSummary	
./Institution/Value	From “Dept. Footer” setting on the Miscellaneous tab of the User Settings.
./Protocol/Value	Name of the stress protocol used. Factory defaults include: <ul style="list-style-type: none"> • BRUCE • Cycle • Modified Bruce • Balke • Ellestad • Naughton • Pharmacological • Low Ramp • Medium Ramp • High Ramp • Astrand • USAF/SAM 2.0 • USAF/SAM 3.3
./Protocol/Label	
./ExamDate/Value	Date of exam, in format displayed to user.
./ExamDateEx/Value	Date of exam, in XML format, yyyy-MM-dd.
./ExamTime/Value	Local time of day when exam began in hh:mm format.
./ExerciseTime/Value	Total exercise time in h:mm:ss format.
./JPoint/Value	Milliseconds from J-point where ST level is measured.
./JPoint/Unit	Always ms (milliseconds).
./JPoint/Label	

XML Tag	Description
./LeadsWith100uV_ST/Value	One Value for every lead having at least 100 uV of ST elevation or depression. Values can be: <ul style="list-style-type: none"> • I • II • III • aVR • aVL • aVF • V1 • V2 • V3 • V4 • V5 • V6
./PVCs/Value	Total number of PVCs detected during the exam.
./PVCs/Label	
./DukeScore/Value	Duke treadmill score when the Bruce exercise protocol is used. Ranges from approximately -57 to 21.
./DukeScore/Label	
./FAI/Value	Functional Aerobic Impairment score, expressed as a percentage. Two values are given separated by a / slash. The first value listed is for a person with a sedentary lifestyle (does not exercise at least once per week sufficient to perspire) and the second value is for an active person (exercises at least once per week sufficient to perspire).
./FAI/Label	
./MaxSpeed/Value	Maximum treadmill speed during exam. Expressed as a number with units (e.g., "5.0 MPH").
./MaxSpeed/Units	<ul style="list-style-type: none"> • MPH = miles per hour • km/h = kilometers per hour
./MaxSpeed/Label	
./MaxGrade/Value	Maximum treadmill grade during exam. Expressed as a number with a percentage sign (e.g., "18.0%").
./MaxGrade/Units	Always %.
./MaxGrade/Label	
./MaxPower/Value	Maximum ergometer power during exam. Expressed as a number.
./MaxPower/Units	Always Watts .
./MaxPower/Label	
./MaxMets/Value	Maximum METs (estimated metabolic equivalents) achieved during test.
./MaxMets/Label	
./MaxHR/Value	Maximum heart rate achieved during exam.
./MaxHR/Units	Always BPM (beats per minute).
./MaxHR/Label	
./MaxSBP/Value	Maximum systolic blood pressure during exam. Expressed as "systolic/diastolic" (e.g. "160/80").
./MaxSBP/SBP	Systolic value.

XML Tag	Description
./MaxSBP/DBP	Diastolic value.
./MaxSBP/Time	Time of measurement, from beginning of the exercise phase. Expressed as h:mm:ss.
./MaxSBP/Units	Always mm Hg (millimeters of mercury).
./MaxSBP/Label	
./MaxDBP/Value	Maximum diastolic blood pressure during exam. Expressed as "systolic/diastolic" (e.g. "160/80").
./MaxDBP/SBP	Systolic value.
./MaxDBP/DBP	Diastolic value.
./MaxDBP/Time	Time of measurement, from beginning of the exercise phase. Expressed as h:mm:ss.
./MaxDBP/Units	Always mm Hg (millimeters of mercury).
./MaxDBP/Label	
./MaxDoubleProduct/Value	Maximum double product (systolic BP * HR) achieved during exam.
./MaxDoubleProduct/Label	
./MaxPercentTargetHR/Value	Maximum percentage of target heart rate achieved during exam.
./MaxPercentTargetHR/Unit	Always %.
./MaxPercentTargetHR/Label	
./MaxST_Elevation/Value	ST level of leading having the most elevation during exam.
./MaxST_Elevation/Units	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./MaxST_Elevation/Lead	Lead having the most ST elevation during the exam.
./MaxST_Elevation/Time	Time elapsed from the beginning of the exam when the maximum ST elevation was measured. Expressed in h:mm:ss format.
./MaxST_Elevation/Label	
./MaxST_Depression/Value	ST level of lead having the most depression during the exam.
./MaxST_Depression/Units	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./MaxST_Depression/Lead	Lead having the most ST depression during the exam.
./MaxST_Depression/Time	Time elapsed from the beginning of the exam when the maximum ST depression was measured. Expressed in h:mm:ss format.
./MaxST_Depression/Label	
./MaxST_ElevationChange/Value	The amount of change measured in lead having the most positive ST change during the exam.
./MaxST_ElevationChange/Units	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./MaxST_ElevationChange/Lead	Lead having the most positive ST change during the exam.
./MaxST_ElevationChange/Time	Time elapsed from the beginning of the exam when the most positive ST change was measured. Expressed in h:mm:ss format.
./MaxST_ElevationChange/Label	
./MaxST_DepressionChange/Value	The amount of change measured in lead having the most negative ST change during the exam.

XML Tag	Description
./MaxST_DepressionChange/Unit	<ul style="list-style-type: none"> • mm = millimeters • uV = microvolts
./MaxST_DepressionChange/Lead	Lead having the most negative ST change during the exam.
./MaxST_DepressionChange/Time	Time elapsed from the beginning of the exam when the most negative ST change was measured. Expressed in h:mm:ss format.
./MaxST_DepressionChange/Label	
./MaxSTHR_Index/Value	The maximum ST/HR index measured during the exam.
./MaxSTHR_Index/Label	
./ReasonsForEnd/Label	
./ReasonsForEnd/Line	<p>One Line per reason. Free text or selected from the customizable pick list configured in ReasonsForEndList.txt. Factory default list is:</p> <ul style="list-style-type: none"> • T1 Target HR • T2 Planned Submax • T3 Chest Pain • T4 Ischemic ECG • T5 Fatigue • T6 Arm Pain • T7 Neck Pain • T8 MD Discretion • T9 Dyspnea Maximum heart rate obtained • T10 Nausea/Headache • T11 Dizziness • T12 Hypotension • T13 Vent Arrhythmia • T14 Atrial Arrhythmia • T15 Claudication • T16 Poor Motivation • T17 Syncope • T18 Completion of Infusion • T19 Completion of Protocol • T20 Other
./Symptoms/Label	
./Symptoms/Line	One Line per symptom. Free text or selected from the customizable pick list configured in SymptomsList.txt .
./Conclusions/Label	
./Conclusions/Line	<p>One Line per line of conclusion block. Free text, template, or selected from the customizable acronym list configured in ConclusionsList.txt. Factory default list is:</p> <ul style="list-style-type: none"> • myo ECG finding suggestive of myocardial ischemia • C1 No ST Changes • C2 Min Depr ST (0.5-0.9 mm) • C3 Mod ST Depr (1.0-1.9 mm) • C4 Marked Depr ST (>=2.0 mm) • C5 ST Elevation During Exercise • C6 Non-Diagnostic TW Abn • C7 Cond Defects Appeared • C8 RBBB Appeared • C9 Inadequate BP Resp (<30)

XML Tag	Description
	<ul style="list-style-type: none"> • C10 Hypertensive BP Response • C11 Inadequate HR c/w Betablocker • C12 Inadequate HR c/w Low level Exercise • C13 Normal Stress Test • C14 Abnormal Stress Test • C15 No Evidence of Ischemia • C16 Chest Pain did not Occur • C17 Atypical Chest Pain • C18 Typical Angina Occured • C19 Exertional Hypotension • C20 Inappropriate Dyspnea • C21 ECG & Sx Typical of CAD • C22 ECG typical of CAD • C23 Sx Typical of CAD • C24 ECG and Sx Suggest CAD • C25 ECG Suggests CAD • C26 Sx Suggest CAD • C27 Markedly abnormal test c/w Extensive CAD • C28 Equivocal Study • C29 Nondiagnostic Test, Inadequate HR Response • C30 Nondiagnostic test, Pt unable to complete • C31 Nondiagnostic Test, Baseline ST depression • C32 Nondiagnostic Test, Conduction Abnormality • C33 Nondiagnostic, Technically Inadequate • C34 Baseline St Abnormality increased during Stress • C35 Baseline ST abnormality unchanged during Stress • C36 Report of Radionuclide study attached • C37 Test within normal limits
./Technician/Value	Name of technician performing the stress exam. Free text or selected from the customizable pick list configured in TechnicianList.txt .
./Technician/Label	
./ReviewingPhysician/Value	Full name of physician reviewing and finalizing the stress report. Free text or selected from the customizable pick list configured in AttendingPhyList.txt .
./ReviewingPhysician/Label	
/StressTest/SummaryTable	
/StressTest/SummaryTable/StageSummaryLine	One StageSummaryLine per stage. One per line of STAGE SUMMARY section of the stress report. Each line reports values at the end of that stage.
./Stage/Value	<p>Name of stage. Values can be:</p> <ul style="list-style-type: none"> • M-LIKAR = Mason-Likar event during pre-exercise • STANDING = Standing event during pre-exercise • HYPERV = Hyperventilation event during pre-exercise • SUPINE = Supine event during pre-exercise • START EXE = end of pre-exercise phase • STAGE 1 = end of stage 1 • STAGE 2 =end of stage 2 • STAGE n = end of stage n, where n is the stage number • PEAK EXE = time and values of peak exercise before entering the recovery phase

XML Tag	Description
	<ul style="list-style-type: none"> RECOVERY = end of 1 minute of recovery. There can be multiple RECOVERY stages, one per minute. END REC = end of recovery
./StageTime/Value	Time when stage ended, expressed as time elapsed from beginning of exercise or recovery phase. Format is: <ul style="list-style-type: none"> PRE-X = pre exercise stage EXE h:mm:ss = exercise stage REC h:mm:ss = recovery stage
./Speed/Value	Treadmill speed at end of stage.
./Speed/Unit	<ul style="list-style-type: none"> MPH = miles per hour km/h = kilometers per hour
./Power/Value	Ergometer workload at end of stage.
./Power/Unit	Always Watts .
./Grade/Value	Treadmill grade at end of stage.
./Grade/Unit	Always % .
./HR/Value	Heart rate at end of stage.
./HR/Unit	Always BPM (beats per minute).
./SystolicBP/Value	Last measured systolic blood pressure in stage.
./SystolicBP/Unit	Always mm Hg (millimeters of mercury).
./DiastolicBP/Value	Last measured diastolic blood pressure in stage.
./DiastolicBP/Unit	Always mm Hg (millimeters of mercury).
./METS/Value	METs (estimated metabolic equivalents) at end of stage.
./DoubleProduct/Value	Double product (systolic BP * HR) at end of stage.
./ST_Level/Lead_I/Value	ST level at end of stage.
./ST_Level/Lead_I/Unit	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./ST_Level/Lead_II/Value	ST level at end of stage.
./ST_Level/Lead_II/Unit	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./ST_Level/Lead_III/Value	ST level at end of stage.
./ST_Level/Lead_III/Unit	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./ST_Level/Lead_aVR/Value	ST level at end of stage.
./ST_Level/Lead_aVR/Unit	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./ST_Level/Lead_-aVR/Value	ST level at end of stage.
./ST_Level/Lead_-aVR/Unit	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./ST_Level/Lead_aVL/Value	ST level at end of stage.
./ST_Level/Lead_aVL/Unit	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./ST_Level/Lead_aVF/Value	ST level at end of stage.
./ST_Level/Lead_aVF/Unit	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./ST_Level/Lead_V1/Value	ST level at end of stage.

XML Tag	Description
./ST_Level/Lead_V1/Unit	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./ST_Level/Lead_V2/Value	ST level at end of stage.
./ST_Level/Lead_V2/Unit	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./ST_Level/Lead_V3/Value	ST level at end of stage.
./ST_Level/Lead_V3/Unit	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./ST_Level/Lead_V4/Value	ST level at end of stage.
./ST_Level/Lead_V4/Unit	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./ST_Level/Lead_V5/Value	ST level at end of stage.
./ST_Level/Lead_V5/Unit	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./ST_Level/Lead_V6/Value	ST level at end of stage.
./ST_Level/Lead_V6/Unit	<ul style="list-style-type: none"> mm = millimeters uV = microvolts
./ST_Slope/Lead_I/Value	Slope of waveform at ST measurement point.
./ST_Slope/Lead_I/Unit	Always mV/s (millivolts per second).
./ST_Slope/Lead_II/Value	Slope of waveform at ST measurement point.
./ST_Slope/Lead_II/Unit	Always mV/s (millivolts per second).
./ST_Slope/Lead_III/Value	Slope of waveform at ST measurement point.
./ST_Slope/Lead_III/Unit	Always mV/s (millivolts per second).
./ST_Slope/Lead_aVR/Value	Slope of waveform at ST measurement point.
./ST_Slope/Lead_aVR/Unit	Always mV/s (millivolts per second).
./ST_Slope/Lead_-aVR/Value	Slope of waveform at ST measurement point.
./ST_Slope/Lead_-aVR/Unit	Always mV/s (millivolts per second).
./ST_Slope/Lead_aVL/Value	Slope of waveform at ST measurement point.
./ST_Slope/Lead_aVL/Unit	Always mV/s (millivolts per second).
./ST_Slope/Lead_aVF/Value	Slope of waveform at ST measurement point.
./ST_Slope/Lead_aVF/Unit	Always mV/s (millivolts per second).
./ST_Slope/Lead_V1/Value	Slope of waveform at ST measurement point.
./ST_Slope/Lead_V1/Unit	Always mV/s (millivolts per second).
./ST_Slope/Lead_V2/Value	Slope of waveform at ST measurement point.
./ST_Slope/Lead_V2/Unit	Always mV/s (millivolts per second).
./ST_Slope/Lead_V3/Value	Slope of waveform at ST measurement point.
./ST_Slope/Lead_V3/Unit	Always mV/s (millivolts per second).
./ST_Slope/Lead_V4/Value	Slope of waveform at ST measurement point.
./ST_Slope/Lead_V4/Unit	Always mV/s (millivolts per second).
./ST_Slope/Lead_V5/Value	Slope of waveform at ST measurement point.
./ST_Slope/Lead_V5/Unit	Always mV/s (millivolts per second).
./ST_Slope/Lead_V6/Value	Slope of waveform at ST measurement point.
./ST_Slope/Lead_V6/Unit	Always mV/s (millivolts per second).

XML Tag	Description
/StressTest/SummaryTable/MinuteSummaryLine	The XML structure of the MinuteSummaryLine is the same as StageSummaryLine described above except for the tags described below. Each MinuteSummaryLine describes the state of the patient at the end of each minute and at manually created events and recorded RPEs.
./Stage/Value	Blank for automatically-created minute lines. Contains the name of the event for manually-created events. If the event is "Shortness of Breath", X-Scribe exports "Short of Breath". For other events, it exports the first 16 characters of the event name. When an RPE is recorded, this is RPE_n , where n is the value on the configured scale of 0-10 or 6-20.

STRESS REQUEST (ORDER) XML

XSD file name: **XScribeInput3_30.xsd**.

File location: on the X-Scribe installation CD, in the **misc** folder.

Stress Request (Order) XML Description

XML Tag	Description
/StressRequest	
./Manufacturer	Leave blank.
./Version	Leave blank.
/StressRequest/PatientDemographics	
./LastName/Value	Patient's last name.
./LastName/Label	Leave blank.
./FirstName/Value	Patient's first name.
./FirstName/Label	Leave blank.
./MiddleName/Value	Patient's middle name.
./MiddleName/Label	Leave blank.
./ID/Value	Patient's primary medical record number.
./ID/Label	Leave blank.
./SecondaryID/Value	Patient's alternate ID. Site-specific use.
./SecondaryID/Label	Leave blank.
./DobEx/Value	Patient's date of birth in XML format, yyyy-MM-dd .
./DobEx/Label	Leave blank.
./Age/Value	Patient's age in years, if DOB is unknown. If DOB is known, age will be calculated at the time of the exam.
./Age/Units	Leave blank.
./Age/Label	Leave blank.
./MaxHR/Value	Leave blank. Calculated by X-Scribe.
./MaxHR/Label	Leave blank.
./TargetHR/Value	Leave blank. Calculated by X-Scribe.
./TargetHR/Label	Leave blank.
./TargetWatts/Value	Target workload for tests with the ergometer.
./TargetWatts/Label	Leave blank.
./Gender/Value	Patient's gender. Values can be: <ul style="list-style-type: none">• Male• Female• Unknown
./Gender/Label	Leave blank.
./Race/Value	Patient's race. Free text. Factory defaults in English are: <ul style="list-style-type: none">• Caucasian• Afro American• Black• Asian• Hispanic• Indian• Other
./Race/Label	Leave blank.
./Height/Value	Patient's height.

XML Tag	Description
./Height/Units	<ul style="list-style-type: none"> • in = inches • cm = centimeters
./Height/Label	Leave blank.
./Weight/Value	Patient's weight.
./Weight/Units	<ul style="list-style-type: none"> • lbs = pounds • kg = kilograms
./Weight/Label	Leave blank.
./Address/Value	Patient's home address. House number and street.
./Address/Label	Leave blank.
./City/Value	Patient's home city.
./City/Label	Leave blank.
./State/Value	Patient's home state.
./State/Label	Leave blank.
./PostalCode/Value	Patient's home postal (zip) code.
./PostalCode/Label	Leave blank.
./Country/Value	Patient's home country.
./Country/Label	Leave blank.
./HomePhone/Value	Patient's home phone number.
./HomePhone/Label	Leave blank.
./WorkPhone/Value	Patient's work phone number.
./WorkPhone/Label	Leave blank.
./ReferringPhysician/Value	Full name of referring physician. Free text.
./ReferringPhysician/Label	Leave blank.
./AttendingPhysician/Value	Full name of the attending physician. Free text.
./AttendingPhysicianLabel	Leave blank.
./Smoker/Value	Values can be: <ul style="list-style-type: none"> • Yes • No • Unknown
./Smoker/Label	Leave blank.
./Diabetic/Value	Values can be: <ul style="list-style-type: none"> • Yes • No • Unknown
./Diabetic/Label	Leave blank.
./HistoryOfMI/Value	Values can be: <ul style="list-style-type: none"> • Yes • No • Unknown
./HistoryOfMI/Label	Leave blank.
./FamilyHistory/Value	Values can be: <ul style="list-style-type: none"> • Yes • No • Unknown
./FamilyHistory/Label	Leave blank.

XML Tag	Description
./PriorCABG/Value	Values can be: <ul style="list-style-type: none"> • Yes • No • Unknown
./PriorCABG/Label	Leave blank.
./PriorCath/Value	Values can be: <ul style="list-style-type: none"> • Yes • No • Unknown
./PriorCath/Label	Leave blank.
./Angina/Value	Values can be: <ul style="list-style-type: none"> • Atypical • Typical • None
./Angina/Label	Leave blank.
./Indications/Line	There is one Line per indication. Lines are free text.
./Indications/Label	Leave blank.
./Medications/Line	There is one Line per medication. Lines are free text.
./Medications/Label	Leave blank.
./Notes/Line	There is one Line per note. Lines are free text.
./Notes/Label	Leave blank.
./MessageID/Value	Values can be used to carry information from the stress order to the stress results. None of the values are displayed to the user on the screen or on the report.
./MessageID/Label	
./OrderNumber/Value	
./OrderNumber/Label	
./BillingCode1/Value	
./BillingCode1/Label	
./BillingCode2/Value	
./BillingCode2/Label	
./BillingCode3/Value	
./BillingCode3/Label	
./ExpansionField1/Value	
./ExpansionField1/Label	
./ExpansionField2/Value	
./ExpansionField2/Label	
./ExpansionField3/Value	
./ExpansionField3/Label	
./ExpansionField4/Value	
./ExpansionField4/Label	

ALTERNATE 1 XML

XSD file name: (not available)

File location: (not available)

NOTE: When ST levels are reported in millimeters, the gain is assumed to be 10 mm/mV.

Alternate 1 XML Description

XML Tag	Description
/Q-Stress_Final_Report	
@UNC	Full path and filename of report in PDF format.
./message_id	
./expansion_field1	
./expansion_field2	
./expansion_field3	
./expansion_field4	
./order_number	
./billing_codes/billing_code	
./machine_id	
./software_version	Product name and version.
/Q-Stress_Final_Report/Summary	
./EvIDProductName	Always X-Scribe .
./EvIDStudyKey	
./EvIDPatientLastName	Patient's last name.
./EvIDPatientFirstName	Patient's first name.
./EvIDPatientMiddleName	Patient's middle name.
./EvIDPatientMRN	Patient's primary medical record number.
./EvIDPatientAccount	
./EvIDPatientSSN	
./EvIDStudyAcqDateISO	Study start date, formatted as yyyy-MM-dd .
./EvIDStudyAcqTimeISO	Local study start time, formatted as HH:mm:ss .
./EvIDStudyDepartment	
./EvIDStudyInstitutionID	
./EvIDStudyInstitution	From "Dept. Footer" setting on the Miscellaneous tab of the User Settings.
./EvIDStudyInstitutionAddress1	
./EvIDStudyInstitutionAddress2	
./EvIDStudyInstitutionCity	
./EvIDStudyInstitutionState	
./EvIDStudyInstitutionZipCode	
./EvIDStudyInstitutionZipCountry	
./EvIDStudySite	
./EvIDStudyAttendingPhysicianEntry	Full name of attending physician.
./EvIDStudyReferringPhysicianEntry	Full name of referring physician.
./EvIDStudyTechnicianEntry	Name of technician.
./EvIDPatientDOBISO	Patient's date of birth, formatted as yyyy-MM-dd .

XML Tag	Description
/EvIDPatientAge	Patient's age in years at time of exam.
/EvIDAgeUnit	Always Years .
/EvIDPatientGender	Patient's gender. Values can be: <ul style="list-style-type: none"> • Male • Female • Unknown
/EvIDPatientHeightValue	Patient's height.
/EvIDPatientHeightUnit	Height units. Values can be: <ul style="list-style-type: none"> • in = inches • cm = kilograms
/EvIDPatientWeightValue	Patient's weight.
/EvIDPatientWeightUnit	Weight units. Values can be: <ul style="list-style-type: none"> • lb = pounds • kg = kilograms
/EvIDPatientAddress1	Patient's home street address.
/EvIDPatientAddress2	
/EvIDPatientCity	Patient's home city.
/EvIDPatientState	Patient's home state.
/EvIDPatientZipCode	Patient's home zip code.
/EvIDPatientCountry	Patient's home country.
/EvIDPatientAddress1Mailing	
/EvIDPatientAddress2Mailing	
/EvIDPatientCityMailing	
/EvIDPatientStateMailing	
/EvIDPatientZipCodeMailing	
/EvIDPatientCountryMailing	
/EvIDPatientAddress1Office	
/EvIDPatientAddress2Office	
/EvIDPatientCityOffice	
/EvIDPatientStateOffice	
/EvIDPatientZipCodeOffice	
/EvIDPatientCountryOffice	
/EvIDPatientPhone	Patient's home phone number.
/EvIDPatientPhoneWork	Patient's work phone number.
/EvIDPatientMedicationEntry	Names of medications. One Entry per medication.
/EvIDFinalRestingHR	Heart rate at end of pre-exercise phase.
/EvIDFinalRestingSysBP	Last measured systolic blood pressure during pre-exercise phase.
/EvIDFinalRestingDiaBP	Last measured diastolic blood pressure during pre-exercise phase.
/EvIDStudyTargetRate	Patient's target heart rate for the exam.
/EvIDStudyMaxPredictedRate	Patient's target heart rate for the exam.
/EvIDFinalPercentMaxHR	Percentage of target heart rate achieved.
/EvIDFinalMaxHR	Maximum heart rate achieved.
/EvIDFinalMaxSysBP	Blood pressure having maximum systolic value, expressed as systolic/diastolic. E.g. "200/110".

XML Tag	Description
./EVIDFinalMaxDiaBP	Blood pressure having maximum diastolic value, expressed as systolic/diastolic. E.g. "180/120".
./EVIDFinalMaxBPStage	
./EVIDProtocol	Name of the stress protocol used. Factory defaults include: <ul style="list-style-type: none"> • BRUCE • Cycle • Modified Bruce • Balke • Ellestad • Naughton • Pharmacological • Low Ramp • Medium Ramp • High Ramp • Astrand • USAF/SAM 2.0 • USAF/SAM 3.3
./EVIDExerciseDevice	Type of exercise device used. Values can be: <ul style="list-style-type: none"> • Treadmill • Ergometer
./EVIDMaxHRxBP	Maximum double product (systolic BP * HR) achieved.
./EVIDFinalOverallWCSlopeValue	
./EVIDFinalOverallWCSlopeLead	
./EVIDFinalOverallWCLevelValue	Worst case ST level (elevation or depression), expressed in millimeters.
./EVIDFinalOverallWCLevelLead	Lead having the worst ST level.
./EVIDFinalTotalExerciseTime	Total exercise time, in h:mm:ss format.
./EVIDFinalMETsAchieved	Maximum METs achieved.
./EVIDLastProtocolStageAchieved	Number of last exercise stage achieved.
./EVIDReasonForTest	Reasons for performing the test. Multiple reasons are separated by spaces.
./EVIDReasonForEndingTest	Reasons for ending the test. Multiple reasons are separated by spaces.
./EVIDTestObservation	Symptoms. Multiple symptoms are separated by spaces.
./EVIDTestConclusion	Test conclusions.
./EVIDExerDevWkldLable	
./EVIDPatientDiagnosisEntry	Diagnosis. One Entry per diagnosis.
./EVIDPatientProcedureEntry	Clinical procedure. One Entry per procedure.
./EVIDPatientRestingECGEntry	
/Q-Stress_Final_Report/Tabular	
/Q-Stress_Final_Report/Tabular/Stage	
@id	Stage. Values can be: <ul style="list-style-type: none"> • REST = pre-exercise • n = exercise stage <i>n</i>, where <i>n</i> is the stage number, starting at 1. • RECOVERY = recovery
@stage_time	Time from beginning of exam, in h:mm:ss format.

XML Tag	Description
./EVIDComment	Name of event. Values can be: <ul style="list-style-type: none"> • Supine • Standing • Chest Pain • Mason-Likar • Shortness of Breath • (free text)
./EVIDExTotalStageTime	Total time in stage, in h:mm:ss format.
./EVIDLogCurrentHR	Heart rate at end of stage.
./EVIDLogCurrentER	
./EVIDLogCurrentBP	Last measured blood pressure in stage, in S/D format. S = systolic, D = diastolic.
./EVIDLogHRxBP	Last measured double product (systolic BP * HR) in stage.
./EVIDExTreadmillSpeed	Treadmill speed at end of stage.
@unit	Units of treadmill speed. Values can be: <ul style="list-style-type: none"> • mph = miles per hour • km/h = kilometers per hour
./EVIDExTreadmillGrade	Treadmill grade at end of stage.
@unit	Always %.
./EVIDSTLevel	ST level at end of stage, in millimeters.
@lead	Lead. Values can be: <ul style="list-style-type: none"> • I • II • III • aVR • aVL • aVF • V1 • V2 • V3 • V4 • V5 • V6
./EVIDSTSlope	ST slope at end of stage, in mV/s.
@lead	Lead. Values can be: <ul style="list-style-type: none"> • I • II • III • aVR • aVL • aVF • V1 • V2 • V3 • V4 • V5 • V6

Below lists possible error codes returned from the DICOM library used by HX-Gate.

DICOM Error Codes

- 1. MC_NORMAL_COMPLETION
- 4000. MC_ALREADY_REGISTERED
- 4001. MC_ASSOCIATION_ABORTED
- 4002. MC_ASSOCIATION_CLOSED
- 4003. MC_ASSOCIATION_REJECTED
- 4004. MC_ATTRIBUTE_HAS_VALUES
- 4005. MC_BUFFER_TOO_SMALL
- 4006. MC_CALLBACK_CANNOT_COMPLY
- 4007. MC_CALLBACK_DATA_SIZE_NEGATIVE
- 4008. MC_CALLBACK_DATA_SIZE_UNEVEN
- 4009. MC_CALLBACK_PARM_ERROR
- 4010. MC_CALLBACK_REGISTERED
- 4011. MC_CANNOT_COMPLY
- 4012. MC_CANT_ACCESS_PROFILE
- 4013. MC_CONFIG_INFO_ERROR
- 4014. MC_CONFIG_INFO_MISSING
- 4015. MC_DDFILE_ERROR
- 4016. MC_DOES_NOT_VALIDATE
- 4017. MC_EMPTY_VALUE
- 4018. MC_END_OF_DATA
- 4019. MC_EXT_INFO_UNAVAILABLE
- 4020. MC_FOUND
- 4021. MC_FUNCTION_UNAVAILABLE
- 4022. MC_INCOMPATIBLE_VR
- 4023. MC_INCOMPATIBLE_VALUE
- 4024. MC_INVALID_APPLICATION_ID
- 4025. MC_INVALID_APPLICATION_TITLE
- 4026. MC_INVALID_ASSOC_ID
- 4027. MC_INVALID_CHARS_IN_VALUE
- 4028. MC_INVALID_COMMAND
- 4029. MC_INVALID_DATA_TYPE
- 4030. MC_END_OF_LIST
- 4031. MC_INVALID_GROUP
- 4032. MC_INVALID_HOST_NAME
- 4033. MC_INVALID_ITEM_ID
- 4034. MC_INVALID_LENGTH_FOR_TITLE
- 4035. MC_INVALID_LENGTH_FOR_VR
- 4036. MC_INVALID_LICENSE
- 4037. MC_INVALID_MESSAG_ID
- 4038. MC_INVALID_MESSAGE_RECEIVED
- 4039. MC_INVALID_PARAMETER_NAME
- 4040. MC_INVALID_PORT_NUMBER
- 4041. MC_INVALID_PRIVATE_CODE
- 4042. MC_INVALID_SERVICE_LIST_NAME
- 4043. MC_INVALID_TAG
- 4044. MC_INVALID_TRANSFER_SYNTAX
- 4045. MC_INVALID_VALUE_FOR_VR
- 4046. MC_INVALID_VALUE_NUMBER
- 4047. MC_INVALID_VR_CODE

4048. MC_LOG_EMPTY
4049. MC_MESSAGE_EMPTY
4050. MC_MESSAGE_VALIDATES
4051. MC_MISSING_CONFIG_PARM
4052. MC_MSGFILE_ERROR
4053. MC_MUST_BE_POSITIVE
4054. MC_NETWORK_SHUT_DOWN
4055. MC_NO_APPLICATIONS_REGISTERED
4056. MC_NO_CALLBACK
4057. MC_NO_CONDITION_FUNCTION
4058. MC_NO_FILE_SYSTEM
4059. MC_NO_INFO_REGISTERED
4060. MC_NO_LICENSE
4061. MC_NO_MERGE_INI
4062. MC_NO_MORE_ATTRIBUTES
4063. MC_NO_MORE_VALUES
4064. MC_NO_PROFILE
4065. MC_NO_REQUEST_PENDING
4066. MC_NON_SERVICE_ATTRIBUTE
4067. MC_NOT_FOUND
4068. MC_NOT_ONE_OF_ENUMERATED_VALUES
4069. MC_NOT_ONE_OF_DEFINED_TERMS
4070. MC_NULL_POINTER_PARM
4071. MC_NULL_VALUE
4072. MC_PROTOCOL_ERROR
4073. MC_REQUIRED_ATTRIBUTE_MISSING
4074. MC_REQUIRED_DATASET_MISSING
4075. MC_REQUIRED_VALUE_MISSING
4076. MC_STATE_VIOLATION
4077. MC_SYSTEM_CALL_INTERRUPTED
4078. MC_SYSTEM_ERROR
4079. MC_TAG_ALREADY_EXISTS
4080. MC_TEMP_FILE_ERROR
4081. MC_TIMEOUT
4082. MC_TOO_FEW_VALUES
4083. MC_TOO_MANY_BLOCKS
4084. MC_TOO_MANY_VALUES
4085. MC_UNABLE_TO_CHECK_CONDITION
4086. MC_UNACCEPTABLE_SERVICE
4087. MC_UNEXPECTED_EOD
4088. MC_UNKNOWN_ITEM
4089. MC_UNKNOWN_SERVICE
4090. MC_VALUE_MAY_NOT_BE_NULL
4091. MC_VALUE_NOT_ALLOWED
4092. MC_VALUE_OUT_OF_RANGE
4093. MC_VALUE_TOO_LARGE
4094. MC_VR_ALREADY_VALID
4095. MC_LIBRARY_ALREADY_INITIALIZED
4096. MC_LIBRARY_NOT_INITIALIZED
4097. MC_INVALID_DIRECTORY_RECORD_OFFSET
4098. MC_INVALID_FILE_ID
4099. MC_INVALID_DICOMDIR_ID

4100. MC_INVALID_ENTITY_ID
4101. MC_INVALID_MRDR_ID
4102. MC_UNABLE_TO_GET_ITEM_ID
4103. MC_INVALID_PAD
4104. MC_ENTITY_ALREADY_EXISTS
4105. MC_INVALID_LOWER_DIR_RECORD
4106. MC_BAD_DIR_RECORD_TYPE
4107. MC_UNKNOWN_HOST_CONNECTED
4108. MC_INACTIVITY_TIMEOUT
4109. MC_INVALID_SOP_CLASS_UID
4110. MC_INVALID_VERSION
4111. MC_OUT_OF_ORDER_TAG
4112. MC_CONNECTION_FAILED
4113. MC_UNKNOWN_HOST_NAME
4114. MC_INVALID_FILE
4115. MC_NEGOTIATION_ABORTED
4116. MC_INVALID_SR_ID
4117. MC_UNABLE_TO_GET_SR_ID
4118. MC_DUPLICATE_NAME
4119. MC_DUPLICATE_SYNTAX
4120. MC_EMPTY_LIST
4121. MC_MISSING_NAME
4122. MC_INVALID_SERVICE_NAME
4123. MC_SERVICE_IN_USE
4124. MC_INVALID_SYNTAX_NAME
4125. MC_SYNTAX_IN_USE
4126. MC_NO_CONTEXT
4127. MC_OFFSET_TABLE_TOO_SHORT
4128. MC_MISSING_DELIMITER
4129. MC_COMPRESSION_FAILURE
4130. MC_END_OF_FRAME
4131. MC_MUST_CONTINUE_BEFORE_READING
4132. MC_COMPRESSOR_REQUIRED
4133. MC_DECOMPRESSOR_REQUIRED
4134. MC_DATA_AVAILABLE
4135. MC_ZLIB_ERROR
4136. MC_NOT_META_SOP
4137. MC_INVALID_ITEM_TRANSFER_SYNTAX
4138. MC_LICENSE_ERROR

