

Doc: 076-8431-00  
Rev: A  
File: 98200-117-IX  
Date: 21 June 2016  
Page 1 of 2

**Equipment Line: Cardiology – Data Management**

**Subject: Sentinel Configuration for CardioExpress SL6A/12A/18A Orders**

**Affected: Sentinel (98200)**

**Applicability:**

Mod Category: N/A  
Urgency: N/A  
Serial Number: N/A  
Time Required: N/A.  
In-Service: N/A

**Dependency:**

Applies only with “98410” series CardioExpress models SL6A, SL12A, SL18A.

**Explanation:**

These instructions describe an additional configuration for the Sentinel Cardiology Information Management System to send a worklist (orders) to “98410” series CardioExpress models SL6A, SL12A, and SL18A. This configuration is in addition to PSN 076-8169-00, CardioExpress Bi-Directional Configuration. It is compatible with the earlier 98400 series CardioExpress SL6/SL12, used in the same installation.

**Implementation**

Before making any changes, make a backup of the current file.

In a Sentinel Network environment, make the change on the Sentinel Server.

In a Sentinel Stand-alone configuration, make the change on the workstation.

**Additional Configuration for Sentinel 9 (all versions):**

| Action  | Expected result  |
|---|--|
| Open <b>testsToQueryResponse.xslt</b> in a text editor in Administrator mode. Default location is “C:\Program Files (x86)\Del Mar Reynolds Medical\Sentinel Messaging Service\Steering Files\testsToQueryResponse.xslt” | Configuration file open  |
| Search for the line containing “<record>”   | Line located   |
| Insert the following text as a new line immediately after the line located above:<br><b>&lt;SerialNO&gt;&lt;xsl:value-of<br/>select="OrderNumber"/&gt;&lt;/SerialNO&gt;</b>   | Line inserted. Note that this line should be present between <record> and the line containing <PatientID>...</PatientID> |
| Save the file, and then restart “Sentinel Messaging Service” from services control panel or restart computer  | Service / computer restarts as necessary   |

**COMPANY CONFIDENTIAL**  
FOR AUTHORIZED USE ONLY  
© Spacelabs Healthcare 2016

**Additional Configuration for Sentinel V10.0.0:**

| Action  | Expected result   |
|---|---|
| Open <b>testsToQueryResponse.xslt</b> in a text editor in Administrator mode. Default location is "C:\Program Files (x86)\Del Mar Reynolds Medical\Sentinel Messaging Service\Steering Files\testsToQueryResponse.xslt" | Configuration file open   |
| Search for the line containing "<record>"   | Line located  |
| Insert the following text as a new line immediately after the line located above:<br><b>&lt;SerialNO&gt;&lt;xsl:value-of<br/>select="OrderID"/&gt;&lt;/SerialNO&gt;</b>   | Line inserted. Note that this line should be present between <record> and the line containing <PatientID> |
| Save the file, and then restart "Sentinel Messaging Service" from services control panel or restart computer  | Service / computer restarts as necessary  |

**Known Issue**

The CardioExpress cannot filter the order list, all orders will appear on every device.

**Parts Required:**

N/A

**Tools and Equipment:**

If necessary, an external disk drive with sufficient capacity to store the backup.

**Parts Deposition:**

N/A

**Reporting:**

N/A

**Attachments:**

N/A

**COMPANY CONFIDENTIAL**  
FOR AUTHORIZED USE ONLY  
© Spacelabs Healthcare 2016